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User's Manual

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Fetus Camera Model M1

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Product package contents list:

Camera	1	Adapter	1
British standard pin	1	European standard pin	1
American standard pin	1	USB charging cable	1
User's Manual	1	Quick Operation Guide	1
Warranty card	1		

Please confirm that the contents of the package are complete. If there is any omission, please contact the customer service staff of Marvoto Company.

CONTENTS >>>

1 Safety Instructions	40
2 Product Introduction	42
3 Product Structure Description	43
4 On/Off and Instructions	44
5 Connecting devices	46
6 Shooting	50
6.1 Professional Mode Shooting	50
6.2 Advanced Mode Shooting	52
6.3 Shooting Class	55
6.4 Photo Browsing, Sharing and Beautifying	65
6.5 Frequency and Duration of Use	66
7 Charging	67
8 Firmware Upgrade	68
9 Storage, Cleaning and Maintenance	70
10 Troubleshooting	72
11 Specifications	74
12 Acoustic Output Parameter Table	75

1 Safety Instructions

This instrument is an internal power supply device, of which the type of electric shock protection is BF type application part.

To avoid possible injury, please read the following safety instructions carefully before using this instrument.

Safety Alert Symbol Description

Symbol	Name	Description
	Prohibition	Behaviour that should be avoided
	Caution	Behaviour that may cause damage to the product or some danger. Please read carefully.
	Attention	General attention for operation and use

 Please use this product in accordance with the laws of the relevant country or region.

 This product is intended for women of pregnancy of 16 weeks and above.

 This product is only used for fetus scanning of the abdomen of pregnant women.

 This product cannot be used for medical diagnosis (e.g. fetus teratology examination).

 Please do not use this product for newborns.

 Do not use if pregnant women have fever and illness.

 This product has requirements on frequency of use and duration of a single time use, so please strictly follow the instructions described in the manual.

 This product is a precision device and should not be disassembled.

- ⚠ Do not break off or press the handle with force to avoid damage.
- ⚠ This product is not completely waterproof. Do not completely immerse the device in water or other liquids.
- ⚠ Be careful not let liquid or other foreign matter into the heat dissipation hole.
- ⚠ Do not allow children to operate or play with the device to avoid accidentally breaking the device or causing injury to children.
- ⚠ Please keep the device clean and avoid oscillation.
- ⚠ Ensure that the product's environment is not subject to strong electromagnetic interference sources such as wireless transmitters, microwave ovens, induction cookers, etc.
- ⚠ Keep the operating environment free from vibration, corrosion, flammable materials, and avoid excessive high or excessive low temperature and humidity.
- ⚠ If this instrument is splashed or has water condenses, stop operation.
- ⚠ Check and ensure that the device has no clear damage before use.
- ⚠ This product must be repaired by an authorized and qualified engineer.
- ⚠ After the product has expired, please return it to the manufacturer for recycling and handling according to local regulations.

2 Product Introduction

The Fetus Camera M1 is a new concept camera designed for pregnant women to take pictures and record dynamic video of fetus.


- ◆ Applicable to: women of pregnancy of 16 weeks and above
- ◆ Applicable parts: abdomen of pregnant women

The Fetus Camera M1 is composed of a device end (hereinafter referred to as "the camera main body") and an APP end (APP name: fetus camera). The camera main body uses the principle of ultrasound imaging. Its core is the Intelligent fetus imaging system, which is responsible for transmitting and receiving ultrasonic signals, signal conversion processing, digital processing and 3D reconstruction processing, etc. The image data output by the camera main body is transmitted to the APP for display via Wi-Fi, and stored, browsed, socially shared, etc., on the APP.

The Fetus Camera has passed the following certifications and tests:

CE ROHS

3 Product Structure Description

1——Connection indicator

2——Mode indicator

3——Power button (with indicator)

4——Ultrasonic sensor

5——Main body

6——Heat dissipation hole

7——L-shaped handle

8——Top cover

9——Type-C USB interface

4 On/Off and Instructions

Indicator Colour and Status Description

Indicator	Colour and Status	Description
Power Indicator	The blue light is always on	The power is sufficient or the power is full at the status of turning on.
	The red light is always on	Charging...
	Slow red flash	Electricity $\leq 20\%$
	Fast red flash	Electricity $\leq 10\%$
	Light off	Turning off or the power is full at the status of turning off.
Mode indicator	The green light is always on	Advance mode/profession mode on
	Yellow and green lights are flashing alternately	Turning on...
	Yellow and green lights are flashing alternately and slowly	Sleeping...
	Light off	Turning off

Connection indicator	Red and blue lights are flashing alternately	Waiting for the connection...
	The blue light is always on	Connection succeeded
	Light off	Turning off

Turning on:

Press and hold the power button for 3 seconds, and release when the mode light "is flashing yellow-green alternately". The camera enters the turning-on program, and when the three lights are all on, the camera is turned on successfully.

Turning off:

Press and hold the power button for 3 seconds, and release when the light is off. The camera will be turned off.

After the power is turned on, the ultrasonic sensor at the bottom of the camera will have an automatic calibration. When the camera is turned on after being laid aside for a long time, the auto calibration time will be prolonged and a short subtle sound may be heard, which is normal.

The camera will automatically turn off in the following situations:

- 1)The camera's power is as low as 5%;
- 2)Improper use causes the camera temperature to reach the warning value of 90 °C.

5 Connecting devices

A "point-to-point" Wi-Fi transmission method is used between the Fetus Camera and the mobile phone.

The mobile phone is directly connected to the camera main body for data communication

After the camera is turned on, it will broadcast a Wi-Fi hotspot. You can start the photo shooting only when the mobile phone connects to such Wi-Fi hotspot and the APP establishes a connection with the device.

◆ Wi-Fi hotspot name: Marvoto_M1_*****

◆ Wi-Fi password: 12345678

Connection procedure

1. Please check the camera and mobile phone power before use to ensure sufficient power.
 2. Please confirm that the mobile phone's system is iOS7.0 or above or Android7.0 or higher.
-

Step 1 Install the "Fetus Camera" APP

Launch the mobile phone browser and scan the QR code below (or the QR code on the product packing box) to download and install the "Fetus Camera"

iOS users can enter the words "Fetus Camera" directly in the APPStore to search and download.

Android users can enter the words "Fetus Camera" in Googleplay and MyApp market to search and download.


Step 2 Registration and login

Launch the APP and go to the registration page. Select the email registration or mobile phone number registration, and follow the page instructions to complete the registration and login.

Step 3 Add a device


After registering and logging in, go to the homepage and click "Add Device" to enter the adding process:

The iOS user shall operate as follows:

Turn ON the unit and Enter "setting"—"Wi-Fi" Open the APP once more. enter the next step after page on the mobile phone. The page will show the red and blue lamp Find and connect to the successful connection. alternately flashing. WIFI named as "Marvoto_ Click "Next". M1_" with the password 12345678.

The Android user shall operate as follows:

Turn on the unit and When the mobile phone Connection succeeded enter the next step after finds the device, select the red and blue lamp the device and click alternately flashing. "Connect Now"

Android users click "Take a Photo" to enter the shooting interface; iOS users click "Start Shooting" button on the homepage to enter the shooting interface.

After the mobile phone is connected to the Fetus Camera, the mobile phone will not be able to access the Internet, so the functions that need to access the Internet cannot be used, such as the "Share" and "Submit" functions in the APP, as well as the Internet functions of other applications on the mobile phone. It is necessary to disconnect from the camera and switch to a normal Wi-Fi connection before you can access the Internet.

6 Shooting

The Fetus Camera has two modes of operation: "Professional Mode" and "Advanced Mode".

The professional mode captures 2D video and the advanced mode captures 3D video.

2D image taken in professional mode – black and white fan-shaped image

3D image taken in advanced mode – real person skin colour image

Photo shooting preparation: Find a comfortable sofa or bed and let the pregnant mother recline on the sofa or bed.

6.1 Professional mode shooting

Find the fetus position on the fan-shaped image

Turn or tilt the camera to find the best shooting angle


Click "Record" to record 2D video

Step 1 Apply couplant

Unscrew the couplant cap and squeeze the couplant evenly onto the pregnant mother's belly. In order to confirm the position of the fetus, the application area can be larger.

Step 2 Place the camera on the belly where the couplant is applied

Holding method


Gently press to make the camera fit better with the abdomen skin, which helps to capture sharper images.


Step 3 Slowly move the camera to find the fetus' position

On the fan-shaped image in black and white, when the position is right, you can see the fetus' image when placing the camera onto the belly.

In order to get a better shooting angle, you need to slowly change the camera angle or move the camera position.

Before rotating

After rotating the camera
90 degrees

Try to adjust the brightness and contrast to get clear fetus image with the best brightness.

Step 4 Click the "Record" button to record the video.


Click the "●" button at the bottom of the interface to start recording 2D video.

Click the "⊙" button to end the recording. Click the "📷" button to capture and save the current frame image to the album.

6.2 Advanced Mode Shooting


In the professional mode, click "Advanced Mode" or slide the menu area to the left to switch to "Advanced Mode".

In the advanced mode, you will first see a black and white fan-shaped image that is the same as the professional mode, because the advanced mode needs to go through "framing" on the fan-shaped image first - find the fetus position, determine the shooting angle, adjust the finder frame, and click "Shooting" to process and generate 3D video.


Step 1 Find the fetus and adjust the "finder frame" to the right position and size

After finding the fetus' position and adjusting the camera angle, you need to adjust the "finder frame":

Red frame - activated and adjustable
White dotted line - cannot be adjusted


White frame - cannot be adjusted
Red dotted line - activated and adjustable

The finder frame changes to "red", indicating that the finder frame is activated in an adjustable state:

- ◆ Move the finder frame: Move the finger in the image area, and the finder frame follows the movement; ;
- ◆ Zoom the finder frame: Two fingers zoom on the screen, and the finder frame follows the zoom.

The dotted line in the middle of the finder frame is a split line, and the image area below the dotted line in the finder frame will be the content of the 3D image.


The image area where the 3D image will be generated and the image outside the frame will not generate a 3D image.

Click the "" icon at the bottom of the page to switch and activate the split line adjustment (continuous click to switch the "finder frame" and "split line" adjustments):


- ◆ Slide up or down the finger on the dotted line, and the dotted line follows it to bend upward or downward;
- ◆ Slide the finger up/down on both ends of the dotted line to move the endpoint position.

Step 2 Click "Shooting" to enter real-time 3D video.


After determining the shooting angle and the position and size of the finder frame, stop moving the camera and maintain the current pressing force. Click the  button at the bottom of the page to process and generate real-time 3D image of the fetus.

Advance Mode preparation

Real-time 3D video


Save the video

3D video image optimization adjustment


The real-time video page provides rotation, light, and filter adjustments. Try adjusting to help you get a better 3D image.

1) Filter adjustment


Try sliding the slider or clicking the "+" and "-" buttons. The image information on the 3D photo will follow the operation to become less or more, and adjust the filter value until you can see the subject clearly.

The filter value is too large and part of the fetus body's data is filtered out.

The filter value is too small, causing too much extra interference information

The filter value is right, so the fetus image is clear, full, and has less interference.

2) Lighting

Click the "💡" icon and there will be a light effect on the 3D image

(similar to "lighting" when shooting in the studio), and you can switch between the lights in different directions by clicking continuously.

Step 3 Click the "End Recording" button to end the recording.

Click the "⊖" button at the bottom of the page to complete the current video recording and enter the "Save" interface, click "↓" to save the video.

6.3 Tips

What is the role of the couplant?

An ultrasonic couplant is a gel made of water as a main component. In the ultrasound examination of the human body, the air between the ultrasonic probe and the skin affects the emission-reception of the ultrasonic wave, therefore the couplant is applied to make the probe and the skin closely fit, and the influence of the air is excluded, which is helpful for collecting clearer images. The couplant also makes it easier for the probe to slide over the skin surface.


After the photo is taken, wipe off the couplant on the belly and the camera with a paper towel.

Can't read the content on the fan-shaped image?



The fan-shaped image has a professional name as "two-dimensional image" or "section image", which is a flat image composed of many light spots of different brightness.

The brighter the area on the fan-shaped image (the whiter the colour), indicating the stronger the ultrasonic echo signal of the human tissue, and the darker the area (the darker the colour), indicating the weaker the ultrasonic echo signal.

In general, the image of the human body's bones, blood vessel walls, fascia, etc., will be brighter on a two-dimensional image. The images of blood, amniotic fluid, urine and other liquids on a two-dimensional image are black, and the muscles, liver and other parenchymal organs or tissues are gray images of varying brightness.

Taking the above picture as an example. This is a two-dimensional image of a fetus of about 26 weeks. The clear black area on the image is amniotic fluid. Because the amniotic fluid is relatively clean and the black and white contrast is clear, the fetus head at the bottom can be easily identified.

Relatively speaking, the fetus' head, hands, feet or legs are relatively easy to distinguish, and the fetus is moving on the two-dimensional image. After seeing several times, the pregnant mother will gradually become familiar with the fetus' head, hands and feet.

How do you shoot a fan-shaped image?

The camera placement angle changes slightly, and the fan-shaped image display content will be different.


Taking the shooting of fetus face as an example, in order to take a three-dimensional front photo of the fetus face, it is necessary to take a longitudinal section image of the fetus face upward (that is, taking the fetus' side face).

 Longitudinal section image

Imagine the fetus lying in front of you and connect the fetus head to the buttocks. The plane image parallel to the line is the fetus' longitudinal section image.

In the picture above, the fetus is head down in the mother's stomach (head facing the mother's leg, the fetus position is the head position), and the face is towards the mother's belly, which is an ideal position to shoot. By placing the camera upright on the mother's belly, you can get the fetus' longitudinal section image.

In actual shooting, the fetus position may not be like the above picture, so you can first find the shape of the fetus head, as shown below:

The white round area with dark interior is the fetus head

Then rotate or tilt the camera angle to quickly find the side face position, as shown below:

If the fetus face is facing the mother's back, the fan-shaped image shows the fetus' back when the camera is placed on the mother's belly. At this time, move the camera to the sides of the mother's body, and a part of the fetus' side face may be photographed. If you can't see the fetus face, it means that the fetus position is not suitable, and you can only shoot again in another time.

💡: How can I take a 3D image of the fetus face?

First of all, the pregnancy weeks should be appropriate.

The best time to take image of the fetus face is pregnancy of 16W-35W. At this time, the fetus head is not too big, while the amniotic fluid in the uterus is sufficient, and the light transmission is good. Thus it is easy to take a clear three-dimensional image of the fetus.

When the pregnancy week is too small (before 16W), for the sake of your fetus' health and safety, try to reduce the shooting times.


When the pregnancy week is too large (after 35W), the fetus has grown up, while the amniotic fluid in the uterus is reduced and the light transmission is weakened. Thus it is difficult to capture the fetus' three-dimensional image.

Second, the fetus position should be appropriate.

The ideal position for shooting is that the fetus face is facing the mother's belly; if the fetus face is facing the mother's back, the fetus position is not suitable, and it is difficult to photograph the fetus face.

Third, "framing" is critical.

The three-dimensional image content is generated by processing a black and white two-dimensional image of continuous multi-frames through a three-dimensional reconstruction algorithm.


Therefore, the content of the three-dimensional image is directly determined by the content presented on the two-dimensional image. If the angle of the two-dimensional image is not good or not clear, the content of the three-dimensional image will not be ideal, so the "framing" of the two-dimensional image is very important.


A good "framing" should meet the following three conditions:

- ◆ The camera angle is correct – you can see the subject in the fan-shaped image;
- ◆ Ample amniotic fluid – black amniotic fluid area around the shooting subject;
- ◆ The finder frame is reasonable – the shooting subject is in the finder frame, and the split line is curved to fit the shooting subject.

Example of shooting a fetus face:

Optimum viewing angle sketch

Three-dimensional image rendering

Schematic diagram of best view angle

Finally, try to adjust the 3D image parameters.

Adjusting 3D image parameters is important for capturing 3D images.

Generally, it can be adjusted as follows:

Step 1 - Rotate the image

Try to rotate the image. Sometimes just rotating the image direction slightly, the angle will be correct, and you can understand the image content.

Step 2 - Turn on the "light" effect

The "light" effect can make the fetus' facial features more solid and help to identify. Try to adjust the light source from different directions until you get the best three-dimension effect of the fetus' facial features on the image.

Step 3 - Adjust the filter

If there is too much information on the image (the flesh colour is full of

image area), try to "turn up" the filter value;

If there is little information on the image (the flesh-coloured image is barely visible), try to "turn down" the filter value;


When shooting the fetus face, the appropriate filter value can make the fetus' facial features clear and full.

Whether you can take the fetus' three-dimensional image, there are many factors such as the size and position of the fetus, more or less amniotic fluid, frequency of fetus movement, etc. The situation of each pregnant mother is different, so if the photo taken is not very good or the fetus position is not suitable for shooting, the pregnant mother can move a little bit or take another time to shoot.

💡 How to quickly find the position of
the fetus' head, hands and feet?

Use fetus position and fetus movement.

- 1) If the pregnant mother has recently gone to the hospital for a check-up, you can initially determine the position of the fetus head according to the fetus position notified by the doctor and try the shooting;
- 2) If you really don't know which part of the stomach to start with, you can find it as follows:


Taking the navel as the centre point, and divide the abdomen into four areas as shown in the figure. Following the order of 1 → 2 → 3 → 4, first look from right to left of the lower abdomen. If the fetus head is not found, then look from the right to the left of the upper abdomen.

- 3) The mother who wants to take a photo of the fetus' hands and feet, just feel the fetus' fetus movement, and put the camera in the position of the fetus movement, then it is very easy to get an image of the fetus' small hands and feet.
-

Advanced mode shooting parameter adjustment technique

Brightness and contrast adjustment

A two-dimensional image with appropriate brightness and obvious contrast between light and dark will help to capture better 3D images.

Brightness value is too high


Brightness value is too low

In general, the default brightness and contrast values of the Fetus Camera can be adapted to the needs of most pregnant mothers.

Finder frame adjustment

The finder frame and split line position must be confirmed before shooting real-time 3D video.

If the finder frame is too small, the 3D image cannot fully present the shooting subject; if the finder frame is too large and too many human tissues are framed, it may interfere with the presentation of the shooting subject. You can just adjust the size of the finder frame to frame the subject.

Split line adjustment

The split line is to exclude the extra interference information above the shooting subject (such as the mother's belly, uterine wall, placenta, etc.), so that the three-dimensional image can show the shooting subject more clearly.

The split line curvature should fit the shooting subject as close as possible, so that the excess information can be removed to the utmost extent.

💡:Common causes of failure of fetus facial 3D image capture

1) The fetus face is close to the placenta: the fetus face is tightly close to the placenta, and there is no gap between the face and the placenta. The split line cannot accurately divide the contour of the face, resulting in a blurred three-dimensional image.

Advice:

- i.Walk around: Move or climb the stairs within the range that you can afford to make the fetus turn and change position;
- ii.Interaction: Pat the belly gently and communicate with the fetus to make he/her change position;
- iii.Adjust parameters: Try adjusting the filter or light effect to see if the image has improved.

2) Insufficient amniotic fluid in the face of the fetus: the distance from the face to the placenta of the fetus or the depth of the amniotic fluid is less than 1 cm.

Because three-dimensional imaging requires the setting of amniotic fluid, if the amniotic fluid is not enough, the facial imaging will be not clear.

Advice:

- i.Supplement amniotic fluid: If the amniotic fluid is not enough in the

early stage, you can supplement it for a while by drinking soup, etc., and you may take an ideal photo later;

ii. Adjust parameters: try to adjust the filter or light effect to see if there is improvement;

iii. Change the shooting angle: Place the camera in a soft position around the fetus head and gently press or tilt the camera angle.

3) The fetus is too active: when the positioning of the finder frame and the split line is completed in the early stage, if the fetus moves or turns during the three-dimensional imaging, it may result in blurred or invisible images;

Advice: The pregnant mother can sit for a while and calm the fetus by appeasement, etc.


4) The fetus is facing away from the outside: If the fetus head back is facing outward (mother's belly), the face is not visible.

Advice:

i. Walk around: Move or climb the stairs within the range that you can afford to make the fetus turn and change position;

ii. Interaction: Pat the belly gently and communicate with the fetus to make him/her change position.

5) There is occlusion in front of the fetus: when there is a hand, foot or umbilical cord in front of the fetus, the face picture taken will be blocked by hands and feet or the umbilical cord. Thus the fetus face will be not full or unclear.

Advice:

i.Walk around: Move or climb the stairs within the range that you can afford to make the fetus turn and change position;

ii.Interaction: Pat the belly gently and communicate with the fetus to make him/her change position.

6) The fetus month is too small: less than 18 weeks and the fetus is not fully developed, so the face is blurred and not 3D enough after imaging; in addition, the internal space of the uterus is sufficient at this stage and the fetus movement is more, so it is difficult to capture the part that you want to shoot.

Advice:

i.Wait for the fetus to mature and then take a face shot to get better results;

ii.Adjust the parameters: Try to adjust the filter or light effect to see if the image is improved;

iii.Change to shoot the fetus' hand or foot. At this time, the fetus is smaller so you can shoot more complete limbs.

7) The fetus month is too large: more than 32 weeks and the fetus head becomes larger. If the position is not ideal, it is not easy to change the position by turning; moreover, most of the pregnant women gradually reduce the amniotic fluid during the post-conceptual age, which may lead to unclear shooting.

Advice:

i.Try to move around or interact with the fetus to see if it can change the position of the fetus;

ii.Adjust the parameters: Try to adjust the filter or light effect to see if the image is improved.

Relatively speaking, the fetus facial shooting requirements are higher. If you can't shoot the ideal three-dimensional photo of the fetus face, you can shoot the fetus' hand, foot or leg, or change to "professional mode".

6.4 Photo Browsing, Sharing and Beautifying

Browse photo

You can browse all saved photos by clicking the "[img]" icon at the bottom of the homepage to switch to the "Album" interface.

Share photo

Click on the photo thumbnail to enter the photo browsing interface. Click the "Share" icon in the upper right corner of the page to share the photo through WeChat, Friends Circle, QQ, Sina Weibo (MicroBlog), Facebook and Twitter.



Photo beautification

The 3D image captured during the shooting can be submitted to the Marvoto cloud server, which will be optimized and processed by Marvoto Cloud IntelligentPhoto processing technology.

Before beautification

After beautification

Users who purchase the camera can enjoy a certain number of free "photo beautification" services. Photo beautification services beyond the agreed range may be charged. Please refer to the information officially published by Marvoto.

6.5 Frequency and Duration of Use

Use frequency:

- ◆ Within the applicable pregnancy weeks (pregnancy of 16 weeks and above), it is recommended to use no more than twice a week;
- ◆ The recommended interval between two usages is 2-3 days.


Single use duration:

- ◆ Do not exceed 10 minutes for each "real use time".
- ◆ Actual use duration: The sum of the actual time that the camera is in contact with the pregnant mother's belly.

7 Charging

- ◆ Install the pin to the adapter:

The product is accompanied with three different pins, so please select the corresponding pin according to the socket standard in your area:

- ◆ Plug one end of the charging cable into the camera's Type-C USB connector, and the other end of the charging cable into the external power outlet.

- ◆ Normally it takes approximately 2.5 hours to complete a single full charge.

Charging precautions:


- ◆ Avoid operating the device during charging;
- ◆ The camera generates heat during charging. Please keep the charging environment ventilated with good heat dissipation;
- ◆ If the camera's charging time is significantly longer, it may be due to poor battery performance or malfunction, please contact our after-sales;
- ◆ Use the original power adapter to charge to avoid accidents when charging.

8 Firmware Upgradew

Check the power of your camera and mobile phone before upgrading to make sure it is fully charged.

Step 1 Download the new firmware version


Click "Firmware Upgrade" on the "Device Details" page to enter the "Firmware Version" page, and click "Download Firmware":

After the firmware download is complete, you need to reconnect the device to start the upgrade.

Step 2 Upgrade the device

After connecting the device, go back to the "Firmware Version" page and click the "Upgrade Device" button to enter the upgrade process:

Please follow the instructions on the page to upgrade.

After the page prompts "Upgrade Successful", please do not perform any operations on the device and wait for the device to turn off itself. After the turning off is completed, it can be turned on again and used normally.

9 Storage, Cleaning and Maintenance

Daily storage and protection

- ◆ The camera is a precision device. When not in use, please put the camera back in the packing box;
- ◆ The sensor at the bottom of the camera is a core component, so please take care not to let it be hit, dropped, etc.
- ◆ Do not scratch/pierce the camera with sharp objects;
- ◆ Avoid contact with corrosive harmful substances;
- ◆ Avoid foreign matter such as water, dust, insects, etc., into the heat dissipation hole;
- ◆ Do not leave the camera in a damp environment;
- ◆ Do not place the camera in an environment where the temperature is above 55°C or below -20°C.

Cleaning and maintenance

Cleaning tools: wiping cloth, detergents, clean water

Cleaning steps:

- ◆ Turn off the camera before cleaning;
- ◆ Clean the dust: Moisten the wiping cloth with water, and wring the water out to gently wipe the surface, and then wipe the water with a dry wiping cloth or paper towel;
- ◆ Clean the dirt: Dip a small amount of detergent with a wiping cloth, and gently wipe the dirt area. After removing the dirt, moisten the wiping cloth with water, and wring the water out to gently wipe again, and then wipe the water with a dry wiping cloth or paper towel;
- ◆ Clean the heat dissipation hole: If dust or foreign matter blocks

the heat dissipation hole, gently pick them out with a cotton swab or toothpick.

Cleaning precautions:

- ◆ Do not spray liquid directly onto the camera to prevent liquid from entering the inside of the camera;
 - ◆ Do not wipe the bottom sensor with alcohol, bleach, ammonium chloride mixture, hydrogen peroxide, etc.
 - ◆ Do not touch solutions or couplants containing mineral oil or lanolin;
 - ◆ If you accidentally spill liquid into the heat dissipation hole, use an air blower to blow dry (with cold air).
-

10 Troubleshooting

Fault	Solution
The three indicators are not all lit up after power on	<p>After normal turning on, the mode indicator will have "yellow and green light flashing alternately". Wait for a few seconds and the connection indicator will light up. If the connection indicator light does not light up after waiting for a long time, please try to turn the device off and restart: press and hold the device "Power" button to tune off the camera and wait for 3 seconds, then press and hold the "Power" button again to see if it is normal; If it is still abnormal, please contact our after-sales for technical support.</p>
No fan-shaped image display after entering the camera interface	<ol style="list-style-type: none"> 1. Click the "Close" button in the upper left corner of the interface to return to the homepage to see if the device is connected. If it shows "Not Connected", please connect the device; 2. If the device is connected, click "Start Shooting" to re-enter the camera interface to see if there is a fan-shaped image display; 3. If the operations described in "1" and "2" are not effective, try to enter the mobile phone system task manager to end the APP operation, then re-launch the APP → connect the device → enter the camera interface to see if there is a fan-shaped image display; 4. If the operations described in "1", "2" and "3" are all not effective, try to turn off the device and restart: press and hold the "Power" button to turn off the device and wait for 3 seconds to restart.

<p>In advanced mode, click on the "" button but real-time 3D video cannot display</p>	<ol style="list-style-type: none"> 1. Please confirm that the device and APP are connected properly. 2. Please click the "Close" button in the upper left corner of the page to exit the camera interface and try to shoot again. 3. If there is still no "3D Video" after trying multiple shots, please turn off the device and restart.
<p>The sensor at the bottom of the camera vibrates too much or makes an abnormal sound</p>	<p>Please stop using the device and contact our after-sales for technical support.</p>
<p>When the camera is plugged in with adapter and powered on, the power light flashes red quickly</p>	<p>Normally, after the adapter is plugged in and the power is on, the device is in the charging state, and the red light of the power indicator is always on; if the power light is "flashing red quickly", indicating "charging abnormality", please launch the APP and connect the device to check the error code, and contact our after-sales for technical Support.</p>

11 Specifications

Product name		Fetus Camera
Product model		M1
Adapter	Input	AC100-240V~, 50/60Hz, 1.2A Max
	Output	DC 5.0V /3000mA or DC 9.0V /2000mA or DC12.0V /2000mA
Battery	Battery type	Rechargeable Li-polymer battery
	Battery capacity	2300mAh
	Battery nominal voltage	7.4V
	Battery charge limit voltage	8.4V
Size and weight	Product size	Length × width × height about 161.5mm×105mm×153.5mm
	Main body net weight (with battery)	About 509g
	Weight with package	About 1235g
WiFi	WiFi band	2.4G, 5G
Use duration (reference)	Continuous working hours	1hours
	Standby duration	4hours
	Charging duration	2.5hours
Mobile phone system support		iOS7.0 above or Android7.0 above systems

12 Acoustic Output Parameter Table

Acoustic Output Parameter Table
(According to the standard IEC 61157–2007)

Parameter System settings Standard used: IEC 61157 Ed2	Mode
	B Focus:50mm scan angle:66° frequency:4.0MHz output power:100%
pr(MPa)	1.798
II _{spta} (mW/cm ²)	11.26
I _{ob} (mW/cm ²)	3.30
Power output(mW)	18.23
Output beam dimensions (Ø) (mm)	39
zp(mm)	36
w12 () (mm) (⊥) (mm)	2.60
	5.88
fawf(MHz)	3.603
pr _r (kHz)	4.098
srr(Hz)	21.93
ztt(mm)	—
zts(mm)	contact
Acoustic output freeze	—
Inclusive modes	—